

852 Titrando



Manual
8.852.8003EN



Metrohm AG

CH-9100 Herisau

Switzerland

Phone +41 71 353 85 85

Fax +41 71 353 89 01

info@metrohm.com

www.metrohm.com

852 Titrand

Manual

Teachware
Metrohm AG
CH-9100 Herisau
teachware@metrohm.com

This documentation is protected by copyright. All rights reserved.

Although all the information given in this documentation has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please send us your comments using the address given above.

Documentation in additional languages can be found on
<http://documents.metrohm.com>.

Table of contents

1	Introduction	1
1.1	The Titrand system	1
1.2	Instrument description	2
1.3	Titration modes – Measuring modes – Dosing commands	3
1.4	About the documentation	4
1.4.1	Symbols and conventions	4
1.5	Safety instructions	5
1.5.1	General notes on safety	5
1.5.2	Electrical safety	5
1.5.3	Working with liquids	6
1.5.4	Flammable solvents and chemicals	7
1.5.5	Recycling and disposal	7
2	Overview of the instrument	8
3	Installation	11
3.1	Setting up the instrument	11
3.1.1	Packaging	11
3.1.2	Checks	11
3.1.3	Location	11
3.2	Connecting a controller	11
3.2.1	Operation	11
3.3	Connecting MSB devices	15
3.3.1	Connecting a dosing device	16
3.3.2	Connecting a stirrer or titration stand	17
3.3.3	Connecting a Remote Box	18
3.4	Connecting USB devices	19
3.4.1	General	19
3.4.2	Connecting a USB hub	20
3.4.3	Connecting a printer	20
3.4.4	Connecting a balance	21
3.4.5	Connecting a PC keyboard (only for operation with Touch Control)	22
3.4.6	Connecting a barcode reader	23
3.5	Setting up the titration vessel	24
3.5.1	General	24
3.5.2	Titration vessel for coulometric KF titration	25
3.5.3	Titration vessel for volumetric KF titration	34
3.6	Connecting sensors	38
3.6.1	Connecting a generator electrode (KFC)	38

3.6.2	Connecting an indicator electrode (KFC)	39
3.6.3	Connecting a polarizable electrode (KFT)	39
3.6.4	Connecting a temperature sensor	40
4	Karl Fischer titration	41
4.1	Coulometric titration	41
4.1.1	Principle of coulometry according to Karl Fischer	41
4.1.2	Working with water standards	41
4.1.3	Sample addition	44
4.1.4	Optimum working conditions	45
4.2	Volumetric titration	47
4.2.1	Principle of the volumetric Karl Fischer titration	47
4.2.2	Endpoint determination	47
4.2.3	Karl Fischer reagents	48
4.2.4	Application of the Karl Fischer titration	48
4.2.5	Working with water standards	48
4.2.6	Sample addition	48
4.2.7	Optimum working conditions	51
5	Operation and maintenance	52
5.1	General notes	52
5.1.1	Care	52
5.1.2	Maintenance by Metrohm Service	52
5.2	Generator electrode	53
5.2.1	Generator electrode without diaphragm	53
5.2.2	Generator electrode with diaphragm	53
5.3	Quality Management and qualification with Metrohm	54
6	Troubleshooting	56
6.1	General	56
6.2	Karl Fischer titration	56
6.2.1		56
6.3	SET titration	58
6.3.1		58
7	Appendix	59
7.1	Remote interface	59
7.1.1	Pin assignment of the remote interface	59
8	Technical specifications	63
8.1	Measuring interface	63
8.1.1	Generator electrode	63
8.1.2	Indicator electrode	63
8.1.3	Temperature	63
8.1.4	Polarizer	64

Index 118

Table of figures

Figure 1	The Titrand system	1
Figure 2	Front 852 Titrand	8
Figure 3	Rear 852 Titrand	9
Figure 4	Connecting the Touch Control	12
Figure 5	Connecting the computer	14
Figure 6	MSB connections	15
Figure 7	Connecting a dosing device	17
Figure 8	Connecting an MSB stirrer	18
Figure 9	Connecting the propeller stirrer to the titration stand	18
Figure 10	Connecting the Remote Box	19
Figure 11	Connecting a printer	21
Figure 12	Schematic configuration of magnetic stirrer, indicator electrode and buret tip (Volumetry) /Generator electrode (Coulometry) during a titra- tion. a) stirring direction clockwise, b) stirring direction counterclock- wise.	24
Figure 13	Mounting the coulometer cell	25
Figure 14	Filling the adsorber tube	26
Figure 15	Equipping the coulometer cell	27
Figure 16	Unscrewing the cover from the indicator electrode	29
Figure 17	Unscrewing the cover from the generator electrode	29
Figure 18	Screwing the electrode cable to the electrodes	30
Figure 19	Mounting the addition and aspiration tube	31
Figure 20	Connecting a generator electrode	38
Figure 21	Connecting an indicator electrode	39
Figure 22	Connecting a polarizable electrode	39
Figure 23	Connecting a temperature sensor	40
Figure 24	Connectors of the Remote Box	59
Figure 25	Pin assignment of remote socket and remote plug	59

1 Introduction

1.1 The Titrand system

The Titrand is the heart of the modular Titrand system. Operation is carried out either by Touch Control with a touch-sensitive screen ("stand-alone titrator") or by a computer with a corresponding software.

A Titrand system can contain numerous kinds of a variety of instruments. The following figure provides an overview of the peripheral devices you can connect to the 852 Titrand.




Figure 1 The Titrand system



Up to three control instruments (Titrandos, Dosing Interface, USB Sample Processor, etc.) can be controlled via USB connection during operation with the 900 Touch Control.

You can request information on special applications in the "Application Bulletins" and "Application Notes", available free of charge through the responsible Metrohm representative. Various monographs on the subjects of titration techniques and electrodes are also available.

Updating the device software is described in the Help for the corresponding PC software.

1.2 Instrument description

The 852 Titrando has the following characteristics:

- **Operation**

Operation is carried out by means of a touch-sensitive Touch Control or with high-performance PC software.

- **MSB connectors**

Four MSB connectors (Metrohm Serial Bus) for connecting dosing devices (Dosimat with exchange unit or Dosino with dosing unit), stirrers, titration stands and Remote Boxes.

- **USB connectors**

Two USB connectors, through which devices such as printers, PC keyboards, barcode readers or additional control instruments (USB Sample Processor, Titrando, Dosing Interface, etc.) can be connected.

- **Measuring interface**

Two measuring interfaces. Measuring interface 1 has one measuring input each for:

- a generator electrode
- a temperature sensor (Pt1000 or NTC)
- a double Pt electrode

Measuring interface 2 has one measuring input each for:

- a temperature sensor (Pt1000 or NTC)
- a polarizable electrode

1.3 Titration modes – Measuring modes – Dosing commands

The 852 Titrande supports the following titration modes, measuring modes and dosing commands:

- **MET**
Monotonic equivalence point titration. The reagent addition is carried out in constant volume steps.
Measuring modes:
 - **Ipol** (voltametric measurement with selectable polarization current)
 - **Upol** (amperometric measurement with selectable polarization voltage)
- **SET**
Endpoint titration at one or two specified endpoints, for volumetric bromine index determination.
 - **Ipol** (voltametric measurement with selectable polarization current)
 - **Upol** (amperometric measurement with selectable polarization voltage)
- **KFT**
Volumetric water content determination according to Karl Fischer.
Measuring modes:
 - **Ipol** (voltametric measurement with selectable polarization current)
 - **Upol** (amperometric measurement with selectable polarization voltage)
- **KFC**
Coulometric water content determination according to Karl Fischer.
Measuring mode:
 - **Ipol** (voltametric measurement with selectable polarization current)
- **BRC**
Coulometric bromine index determination. Determining the amount of double bonds in mineral oils, for example.
Measuring mode:
 - **Ipol** (voltametric measurement with selectable polarization current)



- MEAS

The following measuring modes can be selected for measurements:

- **Ipol** (voltametric measurement with selectable polarization current)
- **Upol** (amperometric measurement with selectable polarization voltage)
- **T** (temperature measurement)

- **Dosing commands**

The following commands for dosing can be selected:

- **PREP** (rinsing the cylinder and tubings of an exchange unit or dosing unit)
- **EMPTY** (emptying the cylinder and tubings of a dosing unit)
- **ADD** (dosing a specified volume)
- **LQH** (carrying out complex dosing tasks with a Dosino)

Parallel determinations by coulometric (KFC or BRC) and volumetric titration (KFT or SET/MET) are not possible with the Touch Control.

1.4 About the documentation



CAUTION

Please read through this documentation carefully before putting the instrument into operation. The documentation contains information and warnings which the user must follow in order to ensure safe operation of the instrument.

1.4.1 Symbols and conventions

The following symbols and formatting may appear in this documentation:

(5-12)	Cross-reference to figure legend The first number refers to the figure number, the second to the instrument part in the figure.
1	Instruction step Carry out these steps in the sequence shown.
Method	Dialog text, parameter in the software
File ► New	Menu or menu item
[Next]	Button or key



WARNING

This symbol draws attention to a possible life-threatening hazard or risk of injury.



WARNING

This symbol draws attention to a possible hazard due to electrical current.



WARNING

This symbol draws attention to a possible hazard due to heat or hot instrument parts.



WARNING

This symbol draws attention to a possible biological hazard.



CAUTION

This symbol draws attention to possible damage to instruments or instrument parts.



NOTE

This symbol highlights additional information and tips.

1.5 Safety instructions

1.5.1 General notes on safety



WARNING

This instrument may only be operated in accordance with the specifications in this documentation.

This instrument has left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.

1.5.2 Electrical safety

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.



WARNING

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.



WARNING

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

Mains voltage



WARNING

An incorrect mains voltage can damage the instrument.

Only operate this instrument with a mains voltage specified for it (see rear panel of the instrument).

Protection against electrostatic charges



WARNING

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Do not fail to pull the mains cable out of the mains connection socket before you set up or disconnect electrical plug connections at the rear of the instrument.

1.5.3 Working with liquids



CAUTION

Periodically check all system connections for leaks. Observe the relevant regulations in respect to working with flammable and/or toxic fluids and their disposal.

1.5.4 Flammable solvents and chemicals




WARNING

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set up the instrument in a well-ventilated location (e.g. fume cupboard).
- Keep all sources of flame far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

1.5.5 Recycling and disposal



This product is covered by European Directive 2002/96/EC, WEEE – Waste from Electrical and Electronic Equipment.

The correct disposal of your old equipment will help to prevent negative effects on the environment and public health.

More details about the disposal of your old equipment can be obtained from your local authorities, from waste disposal companies or from your local dealer.



- | | |
|---|--|
| <p>1 Bottle holder
With holding clamps, for two reagent bottles.</p> | <p>2 "On" LED
Lights up when the Titrando is ready for operation.</p> |
| <p>3 "Gen." LED
Lights up when the Titrando is ready for operation and the generator electrode is connected.</p> | |




Figure 3 Rear 852 Titrand

1 Type plate

Contains specifications concerning supply voltage, instrument type and serial number.

3 Connector (Controller)

For connecting a Touch Control or a PC with installed PC software. Mini DIN, 9-pin.

5 Power socket**2 USB connector (USB 1 and USB 2)**

USB ports (type A) for connecting printer, keyboard, barcode reader, additional Titrando, USB Sample Processor, etc.

4 MSB connector (MSB 1 to MSB 4)

Metrohm Serial Bus. For connecting external dosing devices, stirrers or Remote Boxes. Mini DIN, 9-pin.

6 Measuring interface 1 (Input 1) and measuring interface 2 (Input 2)

Input 1 for coulometry.

Input 2 for volumetry.

3.1 Setting up the instrument

The instrument is supplied in highly protective special packaging together with the separately packed accessories. Keep this packaging, as only this ensures safe transportation of the instrument.

Immediately after receipt, check whether the shipment has arrived complete and without damage by comparing it with the delivery note.

The instrument has been developed for operation indoors and may not be used in explosive environments.

Place the instrument in a location of the laboratory which is suitable for operation, free of vibrations, protected from corrosive atmosphere, and contamination by chemicals.

The instrument should be protected against excessive temperature fluctuations and direct sunlight.

3.2.1 Operation

Two different versions are available for operating the 852 Titrand:

- A Touch Control with touch-sensitive screen. It forms a "stand-alone instrument" together with the 852 Titrand.
- A computer enables operation of the 852 Titrand with the help of a PC software, e.g. *tiamo*.



Take care to ensure that the power supply cable is pulled out of the power socket before either setting up or disconnecting connections between the instruments.

3.2.1.1 Connecting a Touch Control



NOTE

The plug is protected against accidental disconnection of the cable by means of a pull-out protection feature. If you wish to pull out the plug, you will first need to pull back the outer plug sleeve marked with arrows.

Connect the Touch Control as follows:

- 1 ■ Insert the plug of the Touch Control connection cable into the **Controller** socket.




Figure 4 Connecting the Touch Control

- 2 ■ Connect the MSB devices (see Chapter 3.3, page 15).
- Connect the USB devices (see Chapter 3.4, page 19).

- 3** ■ Connect the Titrando to the power supply.

- 4** ■ Switch on the Touch Control.

The Touch Control power supply is supplied through the Titrando. Automatic system tests are performed on both instruments at the time of activation. The **On** LED on the front of the Titrando lights up when the system test has been completed and the instrument is ready for operation.



If you would prefer not to position the Touch Control directly next to the Titrand, then you can lengthen the connection with the 6.2151.010 cable. The maximum connection length permitted is 5 m.

The 852 Titrando requires a USB connection to a computer in order to be able to be controlled by a PC software. Using a 6.2151.000 controller cable, the instrument can be connected directly, either to a USB socket on a computer, to a connected USB hub or to a different Metrohm control device.

Cable connection and driver installation

1 Installing the software

- ## 2 Establishing the cable connections

- The "On" LED on the 852 Titrando is not yet illuminated!

- Connect the instrument to a USB connector (Type A) of your computer (see manual of your computer). The 6.2151.000 cable is used for this purpose.




Figure 5 Connecting the computer

The instrument is recognized. Depending on the version of the Windows operating system used, the driver installation proceeds differently afterwards. Either the necessary driver software is installed automatically or an installation wizard is started.

3 Follow the instructions of the installation wizard.


The "On" LED on the 852 Titrande lights up when the driver installation has been completed and the instrument is ready for operation.

If problems should occur during installation, contact your company's IT support team.



NOTE

The plug on the instrument end of the 6.2151.000 controller cable is protected against accidental disconnection by means of a pull-out protection feature. If you wish to pull out the plug, you will first need to pull back the outer plug sleeve marked with arrows.



Registering and configuring the instrument in the PC software

The instrument must be registered in the configuration of your PC software. Once that has been done, you can then configure it according to your requirements. Proceed as follows:

1 Setting up the instrument

- Start the PC software.
The instrument is automatically recognized. The configuration dialog for the instrument is displayed.
- Make configuration settings for the instrument and its connectors.

More detailed information concerning the configuration of the instrument can be found in the documentation for the respective PC software.

3.3 Connecting MSB devices

In order to connect MSB devices, e.g. stirrers or dosing devices, Metrohm instruments are equipped with up to a maximum of four connectors on what is referred to as the *Metrohm Serial Bus* (MSB). Various kinds of peripheral devices can be connected in sequence (in series, as a "daisy chain") at a single MSB connector (8-pin Mini DIN socket) and controlled simultaneously by the respective control device. In addition to the connection cable, stirrers and the Remote Box are each equipped with their own MSB socket for this purpose.

The following figure provides an overview of the instruments that can be connected to an MSB socket, along with a number of different cabling variations.




Figure 6 MSB connections

The control device determines which peripheral devices are supported.



NOTE

When connecting MSB devices together, the following must be observed:

- Only one device of the same type can be used at a single MSB connector at one time.
- Type 700 Dosino and 685 Dosimat dosing devices cannot be connected together with other MSB instruments on a shared connector. These dosing devices must be connected separately.



CAUTION

Exit the control software before you plug in MSB instruments. When it is switched on, the control device automatically recognizes which device is connected to which MSB connector. The operating unit or the control software enters the connected MSB devices into the system configuration (device manager).

MSB connections can be extended with the 6.2151.010 cable. The maximum connection length permitted is 15 m.

3.3.1 Connecting a dosing device

Four dosing devices can be connected to the instrument (**MSB 1 to MSB 4**).

The types of dosing devices that are supported are:

- 800 Dosino
- 700 Dosino
- 805 Dosimat
- 685 Dosimat

Proceed as follows:

1 Connecting a dosing device

- Exit the control software.
- Connect the connection cable of the dosing device to one of the sockets marked with **MSB** on the rear of the control device.
- Start the control software.




Figure 7 Connecting a dosing device

3.3.2 Connecting a stirrer or titration stand

You can use the following instruments:

- With built-in magnetic stirrer (stirring "from below"):
 - 801 Stirrer
 - 803 Ti Stand
- Without built-in magnetic stirrer (stirring "from above"):
 - 804 Ti Stand with propeller stirrer 802 Stirrer

Connect a stirrer or a titration stand as follows:

1 Connecting the stirrer or titration stand

- Exit the control software.
- Connect the connection cable of the magnetic stirrer or of the titration stand to one of the sockets marked with **MSB** on the rear of the control device.
- 804 Ti Stand only: Connect the propeller stirrer to the stirrer connector (socket with stirrer symbol) of the titration stand.
- Start the control software.




Figure 8 Connecting an MSB stirrer




Figure 9 Connecting the propeller stirrer to the titration stand

3.3.3 Connecting a Remote Box

Instruments that are controlled via remote lines and/or that send control signals via remote lines can be connected via the 6.2148.010 Remote Box. In addition to Metrohm, other instrument manufacturers also use similar connectors that make it possible to connect different instruments together. These interfaces are also frequently given the designations "TTL Logic", "I/O Control" or "Relay Control" and generally have a signal level of 5 volts.

Control signals are understood to be electrical line statuses or electrical pulses (> 200 ms) which display the operating status of an instrument or which trigger or report an event. Sequences on a variety of instruments can thus be coordinated in a single complex automation system. No exchange of data is possible, however.

Proceed as follows:

1 Connecting the Remote Box

- Exit the control software.

-

You can, for example, connect an 849 Level Control (fill level monitoring in a canister) or a 731 Relay Box (switch box for 230/110 volt alternating current sockets and low-voltage direct current outlets). The Remote Box also has an MSB socket at which a further MSB device, e.g. a dosing device or a stirrer, can be connected.

3.4 Connecting USB devices

The 852 Titrando has two USB connectors (type A sockets) for peripheral devices with USB interfaces. The Titrando functions as a USB hub (distributor) no matter how it is operated. If you wish to connect more than two devices to the USB, you can also use an additional, commercially available USB hub.



852 Titrando



3.4.2 Connecting a USB hub

If you wish to connect more than two devices to the USB connector of the 852 Titrando, you can also use an additional commercially available USB hub (distributor). If you operate the 852 Titrando with the help of the Touch Control, then you should use a USB hub with its own power supply.

Connect the USB hub as follows:

- 1 Switch off the Touch Control and/or exit the PC software.
- 2 With the aid of the 6.2151.020 cable, connect the USB connector of the 852 Titrando (type A) with the USB connector of the hub (type B, see manual for the hub).
- 3 Switch on the Touch Control.

The USB hub is recognized automatically.

3.4.3 Connecting a printer

Printers that are connected to the 852 Titrando with Touch Control must meet the following requirements:

- Printer languages: HP-PCL (PCL 3 to 5, PCL 3GUI), Canon BJL Commands or Epson ESC P/2
- Printer resolution: 300 dots/inch or 360 dots/inch (Epson)
- Paper size: A4 or Letter, single-sheet feed.

Connect the printer as follows:

- 1 Switch off the Touch Control.
- 2 With the aid of the 6.2151.020 cable, connect the USB connector of the 852 Titrande (type A) with the USB connector of the printer (type B, see manual for the printer).
- 3 Switch on the printer first, then the Touch Control.
- 4 Configure the printer in the device manager of the Touch Control (see Touch Control manual).



- Operation with a PC software:
 - Connect the balance directly to the serial connector (COM) of the computer. This is usually 9-pin and marked with the symbol **IOIOI**.
- Operation with Touch Control:
 - You will need the 6.2148.050 USB/RS-232 adapter to connect a balance.

Balance	Cable
AND ER, FR, FX with RS-232 interface (OP-03)	6.2125.020 + 6.2125.010
Mettler AB, AG, PR (LC-RS9)	In the scope of delivery for the balance
Mettler AM, PM, PE with interface option 016 or Mettler AJ, PJ with interface option 018	6.2146.020 + 6.2125.010 Also from Mettler: ME 47473 adapter and either ME 42500 hand switch or ME 46278 foot switch
Mettler AT	6.2146.020 + 6.2125.010 Also from Mettler: ME 42500 hand switch or ME 46278 foot switch
Mettler AX, MX, UMX, PG, AB-S, PB-S, XP, XS	6.2134.120

Balance	Cable
Mettler AE with interface option 011 or 012	6.2125.020 + 6.2125.010 Also from Mettler: ME 42500 hand switch or ME 46278 foot switch
Ohaus Voyager, Explorer, Analytical Plus	Cable AS017-09 from Ohaus
Precisa balances with RS-232-C interface	6.2125.080 + 6.2125.010
Sartorius MP8, MC, LA, Genius, Cubis	6.2134.060
Shimadzu BX, BW	6.2125.080 + 6.2125.010

Operation with Touch Control

Connect the balance as follows:

- 1 Plug in the USB plug of the USB/RS-232 adapter at the USB connector of the 852 Titrand.
- 2 Connect the RS-232 interface of the USB/RS-232 adapter with the RS-232 interface of the balance (see table for cable).
- 3 Switch on the Touch Control.
- 4 Switch on the balance.
- 5 Activate the RS-232 interface of the balance if necessary.
- 6 Configure the RS-232 interface of the USB/RS-232 adapter in the device manager of the Touch Control (see Touch Control manual).

3.4.5 Connecting a PC keyboard (only for operation with Touch Control)

The PC keyboard is used as an aid for text and numerical input.

Connect the PC keyboard as follows:

- 1 Insert the USB plug of the keyboard into one of the USB sockets of the 852 Titrando.

- The keyboard is recognized automatically and entered in the device manager.

The barcode reader is used as an aid for text and numerical input. You can connect a barcode reader with USB interface.

Connect the barcode reader as follows:

- The barcode reader is recognized automatically and entered in the device manager.

Program the barcode reader as follows (see also the manual for the barcode reader):

- Specify the desired layout for the keyboard (USA, Germany, France, Spain, German-speaking Switzerland).

This setting must match the setting in the device manager (see the Touch Control manual).

- 4 Program the barcode reader in such a way that the ASCII character 02 (STX or Ctrl B) is sent as the first character. This first character is normally referred to as the "Preamble" or "Prefix Code".
- 5 Program the barcode reader in such a way that the ASCII character 04 (EOT or Ctrl D) is sent as the last character. This last character is normally referred to as the "Postamble", "Record Suffix" or "Postfix Code".
- 6 Exit the programming mode.

3.5 Setting up the titration vessel

3.5.1 General

During the titration, it is important that the solution be well-mixed. The stirring rate should be high enough for a small "vortex" to appear. If the stirring rate is too high, then air bubbles will be aspirated. This results in incorrect measured values. If the stirring rate is too low, then the solution at the electrode will not be correctly mixed. In order to ensure that measurement is carried out in a well-mixed solution following addition of the titrant, the buret tip should be placed in a position where the turbulence is high. In addition, the distance between the addition of the titrant and the electrode should be as large as possible. Also take into account the stirring direction (counterclockwise or clockwise) when positioning electrode and buret tip (see figure below).





Figure 12 Schematic configuration of magnetic stirrer, indicator electrode and buret tip (Volumetry) /Generator electrode (Coulometry) during a titration. a) stirring direction clockwise, b) stirring direction counterclockwise.

1	Magnetic stirrer	2	Indicator electrode
3	Buret tip (Volumetry) / Generator electrode (Coulometry)		

3.5.2.1 Mounting the coulometer cell



Mount the coulometer cell as follows on a titration stand:

- 1 Fix the 6.2047.020 titration vessel holder to the 6.2016.050 support rod.
- 2 Insert the 6.1464.320 titration vessel from above into the titration vessel holder.

3.5.2.2 Coulometer cell – Standard setup

Filling the adsorber tube

Before setting up the coulometer cell the 6.1403.030 adsorber tube has to be filled with 6.2811.000 molecular sieve. Proceed as follows:




Figure 14 Filling the adsorber tube

- 1** Insert a small cotton plug into the bottom of the adsorber tube. Do not pack the cotton too tightly.
- 2** Fill the molecular sieve up to $\frac{3}{4}$ of the height.
- 3** Place a small cotton plug on the molecular sieve. Do not pack the cotton too tightly.
- 4** Seal the adsorber tube with the appropriate cover.



NOTE

Note that the molecular sieve must be replaced at regular intervals. Each time you refill the adsorber tube with molecular sieve, you can, for example, write the date directly on the adsorber tube.

Equipping the coulometer cell




Figure 15 Equipping the coulometer cell

Equip the coulometer cell as follows:

- 1 Place the 6.1903.030 stirring bar in the coulometer cell.
- 2 Cut the 6.2713.0x0 ground-joint sleeves to the correct length and attach them to the ground joints of the inserts (electrodes, adsorber tube, etc.).

Take care to ensure that the edges of the ground-joint sleeves are cut to size cleanly and that there are no fringes. The ground-joint sleeves are not permitted to protrude at the lower edge of the ground-joint opening.

- 3** Insert the 6.1403.030 adsorber tube into the generator electrode.



- 4** Insert the 6.0345.100 generator electrode without diaphragm or the 6.0344.100 generator electrode with diaphragm together with the adsorber tube into the large ground-joint opening at the rear.
- 5** Insert the 6.0341.100 indicator electrode into the left ground-joint opening.
- 6** Place the 6.1448.020 septum on the front opening of the coulometer cell and screw it shut with the 6.2701.040 screw cap.

Tighten the screw cap only enough so that it seals. The septum is not permitted to bend.

Filling the coulometer cell (generator electrode with diaphragm)

Proceed as follows when using a generator electrode with a diaphragm:

- 1** Fill approximately 5 mL of catholyte into the generator electrode.
- 2** Fill approximately 100 mL of anolyte into the coulometer cell with the aid of the 6.2738.000 funnel. The level of the anolyte should be roughly 1 - 2 mm above the level of the catholyte.
- 3** Close the remaining ground-joint opening on the right with the 6.1437.000 ground-joint stopper (with ground-joint sleeve attached).

Filling the coulometer cell (generator electrode without diaphragm)

Proceed as follows when using a generator electrode without a diaphragm:

- 1 Fill approximately 100 mL of reagent into the coulometer cell with the aid of the 6.2738.000 funnel.
- 2 Close the remaining ground-joint opening on the right with the 6.1437.000 ground-joint stopper (with ground-joint sleeve attached).

Screwing on the electrode cable

Screw on the electrode cable as follows:

- 1** Unscrew the cover of the indicator electrode.




Figure 16 Unscrewing the cover from the indicator electrode

- 2** Unscrew the cover of the generator electrode.




Figure 17 Unscrewing the cover from the generator electrode

- 3** Tighten the 6.2104.020 electrode cable to the indicator electrode.

- 4** Tighten the 6.2104.120 electrode cable to the generator electrode.



Mark the screw head of the electrode cable. This prevents you from mixing up the indicator and generator electrode.

3.5.2.3 Coulometer cell with addition and aspiration tube (utilization with Ti Stand)





Figure 19 Mounting the addition and aspiration tube

Insert the addition and aspiration tube as follows into the coulometer cell:

- 1** Attach the 6.2713.000 ground-joint sleeve that has been cut to size to the ground joint of the 6.1437.000 stopper.
- 2** Insert the stopper into the 6.1439.010 addition and aspiration tube.
- 3** Attach the 6.2713.000 ground-joint sleeve that has been cut to size to the ground joint of the addition and aspiration tube.
- 4** Insert this assembly into the ground-joint opening.
- 5** Connect the tubing for the reagent addition at the upper connector of the addition and aspiration tube (5).
- 6** Connect the tubing for the aspiration of the coulometer cell at the lower connector of the addition and aspiration tube (6).




When samples release their water only slowly or only at higher temperatures, the oven method is used. The sample is heated in a KF oven (e.g. *860 KF Thermoprep*) and the water that is released is transferred to the coulometer cell with a carrier gas. A detailed description of setting up the coulometer cell with the KF oven can be found in the respective manual.

If a large number of samples have to be processed, the determination of the water content can be automated with the aid of a sample changer with oven module (e.g. *874 Oven Sample Processor*). A detailed description of setting up the coulometer cell with the sample changer can be found in the respective manual.

3.5.3 Titration vessel for volumetric KF titration

Mounting the KF titration cell



Proceed as follows:


- 1 Screw the 6.2013.010 clamping ring tightly to the support rod.
- 2 Fix the 6.1414.030 vessel lid of the KF titration cell (with correctly inserted sealing ring from the 6.1244.040 sealing set) to the support rod. Keep the locking lever pressed down until it can be released at the desired position.
- 3 Fasten the 6.1415.220 (or 6.1415.250) titration vessel with a 6.1903.020 (or 6.1903.030) stirring bar inside on the vessel lid. Fold back the holding bracket upwards while doing so. The markings on the vessel lid and on the plastic ring must be aligned above one another. Afterwards, press the holding bracket downwards in order to fix the titration vessel. The levers of the holding bracket must enclose the pins of the plastic ring on the titration vessel in order to ensure a secure hold.



- Once the height of the KF titration cell has been adjusted correctly, the entire cell can be raised and swiveled as required by pressing the locking lever.


Filling the adsorber tube

Before insertion, the 6.1403.040 adsorber tube must be filled with the 6.2811.000 molecular sieve. Proceed as follows:



- 2** Fill the molecular sieve up to the $\frac{3}{4}$ level.

Inserting the dosing tip, aspiration tip and buret tip



Proceed as follows:


- 1 Place the three screw nipples from 6.2730.030 (including O-rings, but without stoppers) in the rear openings of the vessel lid.
- 2 Insert the 6.1543.110 dosing tip through the screw nipple in the middle rear opening.
- 3 Insert the 6.1543.120 aspiration tip through the screw nipple in the right rear opening.

When solvent is aspirated, the end of the aspiration tip must touch the vessel base, but it must not inhibit the action of the stirring bar.

The aspiration tip can, if needed, be pulled out of the solvent.
- 4 Insert the 6.1543.060 buret tip through the screw nipple in the left rear opening.
- 5 Screw the 6.1805.200 PTFE M8 tubing of the aspiration bottle onto the aspiration tip.

- 6 Screw the 6.1805.200 PTFE M8 tubing of the solvent bottle onto the dosing tip.

Inserting electrode, adsorber tube and septum stoppers



Proceed as follows:

- 1 Introduce the 6.2730.020 septum stopper (with septum inserted) into the front opening of the vessel lid.
- 2 Insert the O-rings of the electrode and of the adsorber tube into the middle openings of the vessel lid.
- 3 Screw the two 6.2730.010 screw nipples into the openings with the O-rings. Do not screw too tightly.
- 4 Introduce the 6.0338.100 double Pt electrode into the left-hand opening and then tighten the screw nipple until it seals.
- 5 Screw the 6.2104.020 electrode cable tightly onto the electrode.

3.6.2 Connecting an indicator electrode (KFC)

Connect the indicator electrode as follows:

- 1** Plug the electrode plug into the **Ind.** socket of the 852 Titrand.




Figure 21 Connecting an indicator electrode



NOTE

The electrode cable is protected against accidental disconnection of the cable by means of a pull-out protection. If you wish to pull out the plug again, you will first need to pull back the outer plug sleeve.

3.6.3 Connecting a polarizable electrode (KFT)

Connect the polarizable electrode as follows:

- 1** Plug the electrode plug into the **Pol.** socket of the 852 Titrand.




Figure 22 Connecting a polarizable electrode



NOTE

The electrode cable is protected against accidental disconnection of the cable by means of a pull-out protection. If you wish to pull out the plug again, you will first need to pull back the outer plug sleeve.

3.6.4 Connecting a temperature sensor

A temperature sensor of the Pt1000 or NTC type can be connected to the **Temp.** connector.

Connect the temperature sensor as follows:

- 1** Insert the plugs of the temperature sensor into the **Temp.** sockets of the 852 Titrand.




Figure 23 Connecting a temperature sensor



NOTE


Always insert the red plug into the red socket. This is the only way that shielding against electrical interference can be ensured.

4 Karl Fischer titration

4.1 Coulometric titration

4.1.1 Principle of coulometry according to Karl Fischer

The **coulometric Karl Fischer titration** is a variation of the classic water content determination method according to Karl Fischer. The conventional method works with a methanolic solution of iodine, sulfur dioxide and a base as buffer substance. If an aqueous sample is titrated, then several reactions take place that can be summarized in the following sum equation:



According to the equation above the I_2 reacts quantitatively with H_2O . This chemical equation serves as a basis for the water content determination.

With the **coulometric Karl Fischer titration**, the necessary iodine is directly and electrochemically generated in the electrolyte containing iodine ("electronic buret"). Between the amount of electric charge and the amount of generated iodine, there is a strictly quantitative relationship, which is used for high-precision dosing of the iodine. Because the coulometric Karl Fischer method is an **absolute determination**, no titer needs to be determined. It must only be ensured that the reaction generating the iodine runs with a 100% current efficiency. All of the reagents available today ensure this.

The endpoint indication is effected voltametrically by modulating an alternating current of constant strength to a double Pt electrode. This results in a voltage differential between the Pt wires. This is drastically reduced as soon as even the slightest amounts of free iodine are present. This circumstance is used for detecting the endpoint of the titration.

4.1.2 Working with water standards

4.1.2.1 Certified water standards

Commercially available, certified water standards with water contents of $1.00 \pm 0.003 \text{ mg/g}$ and/or $0.10 \pm 0.005 \text{ mg/g}$ should be used for validating the instrument as a whole, integrated system.



NOTE

The 1.0 mg/g water standard is easier to handle and therefore to be preferred.

Table 1 Recommended sample size ranges

Water standard 1.0 mg/g	0.2 - 2.0 g
Water standard 0.1 mg/g	0.5 - 5.0 g

4.1.2.2 Practical recommendations

For validation, it is essential to work very accurately. In order to minimize any measurement inaccuracies that could occur, the sample preparation and the sample processing should proceed in accordance with a defined scheme:

- 1 Put on gloves (as always with Karl Fischer titration).
- 2 Use a clean syringe.



NOTE

If you are working with the 0.1 mg/g water standard, then you must use a glass syringe. If you are working with the 1.0 mg/g water standard, then you may use either a plastic syringe or a glass syringe.

- 3 Take a new ampoule of water standard and shake it briefly.
- 4 With a folded paper towel held between thumb and index finger, break open the ampoule at the marking.
- 5 Draw approx. 1 mL of the water standard into the syringe.
- 6 Pull the piston of the syringe up to the end and swing the syringe back and forth somewhat.

The inside of the syringe is being rinsed by water standard and freed from water contamination.
- 7 Dispose of the used water standard in a waste bottle.

- 8** Draw the rest of the water standard into the syringe, aspirating as little air as possible.
- 9** Push out any air bubbles that may be present in the syringe.
- 10** Wipe off the needle with a lint-free paper towel and cover it with the appropriate cap.
- 11** Place the syringe on the balance and press **[TARA]**.
- 12** As soon as the drift on the 852 Titrando is stable, take the syringe by hand, press **[START]** and inject approx. 1 mL of the water standard through the septum.

There are two possibilities:

 - Version 1:
Inject the water standard without immersing the needle in the reagent liquid. If a little drop remains on the end of the needle, it must be aspirated back before pulling the needle out of the septum.
The water standard should not be sprayed from the syringe onto the electrode nor onto the wall of the coulometer cell.
 - Version 2:
Inject the water standard directly under the surface of the reagent liquid.
Take care to ensure that you do not aspirate any liquid when you withdraw the syringe from the reagent liquid.
- 13** Close the syringe with the same cap and place it back on the balance.
- 14** If you have connected a balance to the Titrando, you may transmit the sample size directly from the balance.
- 15** The next determination can be started as soon as the determination has been finished and the coulometer cell has been conditioned (drift stable) again.

Highly volatile samples and samples of low viscosity should be cooled before sampling. Doing so avoids losses while working. The syringe must, however, not be cooled directly, as condensation could form. For the same reason, no air may be aspirated into a syringe into which a cooled sample has been aspirated beforehand.

Samples of high viscosity can be thinned by heating. The syringe must be heated as well. The same effect can also be achieved by diluting with a suitable solvent. In this case, the water content of the solvent has to be determined and subtracted as blank value.

Pastes and fats can be added to the coulometer cell with a syringe without needle. You can use the ground-joint opening for this. If you also wish to aspirate, you can use the opening with the septum stopper. The best way for you to determine the sample amount is to reweigh the sample.

If samples contain only **traces of water**, then the syringe has to be pre-dried well. If possible, the syringe should be rinsed with the sample solution by filling in and discarding solution several times.

4.1.3.3 Working with solid samples

If possible, solid samples are to be extracted or dissolved in a suitable solvent. The resulting solution is injected, and a blank value correction for the solvent has to be made.

If no suitable solvent can be found for a solid sample, or if the sample reacts with the Karl Fischer reagent, then a Karl Fischer oven should be used.

If solid samples are added directly into the coulometer cell, then the generator electrode without diaphragm should be used. The samples can be added through the ground-joint opening or through the side opening. While doing so, take care to ensure that

- the sample releases its moisture completely.
- no side reaction with the Karl Fischer reagent takes place.
- the surfaces of the electrodes are not covered by the sample substance (incomplete KF reaction).
- the Pt grid of the generator electrode does not become damaged.
- the Pt wires of the indicator electrode do not become damaged.

4.1.4 Optimum working conditions

4.1.4.1 General

When a coulometer cell that has been well dried-out beforehand is put into operation with a generator electrode without diaphragm, the basic drift will be reached within approx. 30 minutes. It is recommended that the coulometer cell be repeatedly and carefully shaken during this time.

For generator electrodes with diaphragm, you should expect a preparation time of approx. 2 hours.

To obtain precise determinations of amounts of water smaller than 100 μg , it may also be of advantage to condition the coulometer cell overnight before using it.

4.1.4.2 Drift

A constant drift in the range of $\leq 4 \mu\text{g}/\text{min}$ is all right. Lower values, however, are quite possible. Higher but stable values will still produce good results, because the drift can be compensated.

A constantly high drift can be caused by water-containing deposits in inaccessible parts of the coulometer cell. In these cases, shaking the cell can reduce the value. Ensure that there are no drops above the liquid level in the coulometer cell.

If you are working with a generator electrode with diaphragm, shake the cell only hard enough that the catholyte and anolyte do not mix with one another. If the drift remains too high for a prolonged time, even after shaking the cell, then the electrolyte solutions should be replaced. The catholyte should be replaced once per week.

A wet catholyte can be another reason for the excessively high drift. The wet catholyte can be dried with a KF single-component reagent.

When you work with a Karl Fischer oven, a drift \leq of 10 $\mu\text{g}/\text{min}$ is all right. The drift depends on the gas flow (the smaller the gas flow, the lower the drift) and on the humidity of the surroundings.

4.1.4.3 Reagent replacement

The electrolyte solutions must be replaced in the following cases:

- The coulometer cell is too full.
- The KF reagent has reached its capacity limit.
- The drift is too high, and cannot be reduced by shaking the coulometer cell.
- A two-phase-mixture is being formed in the coulometer cell; in this case it is also possible to aspirate the sample phase only.

Exhausted electrolyte solution is best disposed of by aspiration. To do this, you can use, for example, an *803 Ti Stand* with built-in membrane pump. An advantage is, that the coulometer cell does not have to be disassembled.

In the event of severe contamination, the coulometer cell can be rinsed with a suitable solvent which is also aspirated.

In the case of the generator electrode with diaphragm, the catholyte should be replaced once per week. Longer use can cause blackening and

yellow precipitates in the cathode chamber. A discernable stench is also a sign of having used the catholyte for too long.

4.1.4.4 Indicator electrode

A contaminated indicator electrode can be carefully cleaned with an abrasive agent (6.2802.000 polishing set or toothpaste). After the cleaning, rinse with ethanol.

The two Pt wires of the indicator electrode should run as parallel as possible to one another. Check the Pt wires before inserting the electrode.

4.2 Volumetric titration

4.2.1 Principle of the volumetric Karl Fischer titration

The **volumetric Karl Fischer titration** is the classic method of water content determination. It works with a methanolic solution of iodine, sulfur dioxide and a base as buffer substance. If an aqueous sample is titrated, then several reactions take place that can be summarized in the following sum equation:



According to the equation above, I_2 reacts quantitatively with H_2O . This chemical equation serves as a basis for the water content determination.

The classic Karl Fischer method has undergone continuous further development in recent years. This development did not concern itself solely with the refining and automation of reagent dosing, but also with improving endpoint indication and reagents. The disadvantage of this method is that the reagents are not completely stable. This means that the titer has to be re-determined on a regular basis.

4.2.2 Endpoint determination

The titration endpoint is determined by an electrometric indication method. Double Pt wire electrodes or double Pt ring electrodes are used for this purpose. A distinction is made between the following two indication methods:

Biamperometric indication (Upol)

A constant potential is applied to the electrodes and the resulting current is measured.

Bivoltametric indication (Ipol)

A direct or alternating current is applied between the two electrodes and the resulting potential is measured.

Water content of the sample	KF reagent 1	KF reagent 2	KF reagent 5
1.0%	0.05 - 0.45	0.1 - 0.9	0.25 - 2.25
5.0%		0.02 - 0.18	0.05 - 0.45
10.0%			0.03 - 0.22
25.0%			
50.0%			

Table 4 Approximate sample size in grams (10 mL buret)

Water content of the sample	KF reagent 1	KF reagent 2	KF reagent 5
0.5%	0.2 - 1.8	0.4 - 3.6	
1.0%	0.1 - 0.9	0.2 - 1.8	0.5 - 4.5
5.0%	0.02 - 0.18	0.04 - 0.36	0.1 - 0.9
10.0%		0.02 - 0.18	0.05 - 0.45
25.0%			0.02 - 0.18
50.0%			0.02 - 0.09

Table 5 Approximate sample size in grams (20 mL buret)

Water content of the sample	KF reagent 1	KF reagent 2	KF reagent 5
0.5%	0.4 - 3.6		
1.0%	0.2 - 1.8	0.4 - 3.6	
5.0%	0.04 - 0.36	0.08 - 0.72	0.2 - 1.8
10.0%	0.02 - 0.18	0.04 - 0.36	0.1 - 0.9
25.0%		0.02 - 0.14	0.04 - 0.36
50.0%			0.02 - 0.18

KF reagent 1: 1 mL KF reagent reacts with around 1 mg H₂O

KF reagent 2: 1 mL KF reagent reacts with around 2 mg H₂O

KF reagent 5: 1 mL KF reagent reacts with around 5 mg H₂O

4.2.6.2 Working with liquid samples

Liquid samples are added with a syringe. The samples can be injected two different ways:

- One uses a syringe with a long needle, which one immerses in the reagent during the injection.
- One uses a syringe with a short needle and aspirates the last drops back into the needle.

The best way for you to determine the injected sample amount is to reweigh the sample.

Glass syringes should be used for the **determination of traces and validations**. We recommend obtaining these from a specialized syringe manufacturer.

Highly volatile samples and samples of low viscosity should be cooled before sampling. Doing so avoids losses while working. The syringe must, however, not be cooled directly, as condensation could form. For the same reason, no air may be aspirated into a syringe into which a cooled sample has been aspirated beforehand.

Samples of high viscosity can be thinned by heating. The syringe must be heated as well. The same target can be reached by diluting with suitable solvents. In this case, the water content of the solvent has to be determined and subtracted as blank value.

Viscous samples can be added to the measuring cell with a syringe without needle. You can use the ground-joint opening for this. The best way for you to determine the added sample amount is by reweighing the sample.

4.2.6.3 Working with solid samples

If possible, solid samples are to be extracted or dissolved in a suitable solvent. The resulting solution is injected, during which a blank value correction for the solvent must be carried out.

If no suitable solvent can be found for a solid sample, or if the sample reacts with the Karl Fischer reagent, then a Karl Fischer oven should be used.

If solid samples have to be directly added to the titration cell, they can be inserted through the ground-joint opening. While doing so, take care to ensure that

- the sample releases its moisture completely.
- no side reaction with the Karl Fischer reagent takes place.
- the surface of the electrode is not covered by the sample substance (incomplete KF reaction!).
- the Pt wires of the indicator electrode do not become damaged.

4.2.7 Optimum working conditions

4.2.7.1 Drift

A constant drift in the range of $\leq 10 \mu\text{L}/\text{min}$ is all right. Lower values are, however, quite possible. Higher but stable values will still produce good results, because the drift can be compensated.

A constantly high drift can be caused by water-containing deposits in inaccessible parts of the titration vessel. In these cases, shaking the titration vessel can reduce the value. Make sure that there are no drops above the liquid level in the titration vessel.

When you work with a Karl Fischer oven, a drift $\leq 10 \mu\text{L}/\text{min}$ is all right. The drift depends on the gas flow (the smaller the gas flow, the lower the drift).

4.2.7.2 Reagent replacement

The electrolyte solution must be replaced in the following cases:

- The titration vessel is too full.
- The drift is too high, and cannot be reduced by shaking the titration vessel.

Exhausted electrolyte solution is best disposed of by aspiration. To do this, you can use, for example, an *803 Ti Stand* with built-in membrane pump. An advantage is that the titration vessel does not have to be disassembled.

In the event of severe contamination, the titration vessel can be rinsed with a suitable solvent which is also aspirated.

4.2.7.3 Indicator electrode

A new indicator electrode can take a certain warm-up time for forming the surface. During this time unexpectedly long titration times and high measurement results can occur. This phenomenon will, however, disappear after a short time of use. In order to accelerate the setting of a new indicator electrode, the 852 Titrand can be conditioned e.g. over night.

A contaminated indicator electrode can be carefully cleaned with an abrasive agent (6.2802.000 polishing set or toothpaste). After the cleaning, rinse with ethanol.

The two Pt wires of the indicator electrode should run as parallel as possible to one another. Check the Pt wires before inserting the electrode.

5 Operation and maintenance

5.1 General notes

5.1.1 Care

The 852 Titrando requires appropriate care. Excess contamination of the instrument may result in functional disruptions and a reduction in the life-time of the otherwise sturdy mechanics and electronics.

Spilled chemicals and solvents should be removed immediately. Above all, the plug connections on the rear of the instrument (in particular the power socket) should be protected from contamination.



CAUTION

Although this is largely prevented by design measures, the power plug should be unplugged immediately if aggressive media have found their way into the interior of the instrument to prevent serious damage to the instrument electronics. In such cases, Metrohm Service must be informed.

5.1.2 Maintenance by Metrohm Service

Maintenance of the 852 Titrando is best carried out as part of annual service, which is performed by specialist personnel from Metrohm. If you are frequently working with caustic and corrosive chemicals, a shorter maintenance interval could be necessary.

Metrohm Service offers every form of technical advice for maintenance and service of all Metrohm instruments.

Metrohm representative regarding the precise terms and conditions involved in concluding a corresponding maintenance agreement. Further information on this can be found in the brochure "**Metrohm Care Contracts – Protect your investment the smart way**" available from your local Metrohm representative.

6 Troubleshooting

6.1 General

Problem	Cause	Remedy
The "On" LED is not illuminated, even though the instrument is connected to the power supply.	<i>The Touch Control or the computer has not been switched on yet or the plugs are not correctly plugged in.</i>	<ol style="list-style-type: none"> 1. Check the plug connections. 2. Switch on the Touch Control or the computer.

6.2 Karl Fischer titration

Problem	Cause	Remedy
The drift is very high during conditioning.	<i>The titration cell is leaking.</i>	<ul style="list-style-type: none"> Check the seals and the septum. Replace if necessary. Replace the molecular sieve.
The drift becomes greater after each titration.	<i>The sample releases water very slowly.</i>	<ul style="list-style-type: none"> Adjust the method. Add solubility promoter. Increase the temperature (possibly using a KF oven). See technical literature.
	<i>A side reaction is taking place.</i>	<ul style="list-style-type: none"> Use special reagents. Adjust the method (increase/decrease the temperature, external extraction). See technical literature.
	<i>The pH value is no longer in the optimum range.</i>	Add buffer (see technical literature).
The titration will not be finished.	<i>The titration cell is leaking.</i>	<ul style="list-style-type: none"> Check the seals and the septum. Replace if necessary. Replace the molecular sieve.
	<i>The minimum increment is too low.</i>	Select the user-defined titration rate and increase the minimum volume increment (see manual/help of the software used).
	<i>The stop criterion is unsuitable.</i>	Adjust the control parameters (see manual/help of the software used): <ul style="list-style-type: none"> Increase the stop drift.

Problem	Cause	Remedy
		<ul style="list-style-type: none"> Select a short delay time.
	<i>See also: The drift becomes greater after each titration.</i>	
The sample is over-titrated.	<i>The increments at the end of the titration are too high.</i>	<ul style="list-style-type: none"> Select the user-defined titration rate and reduce the dosing rate (see manual/help of the software used). The following experiment provides a clue for the optimum dosing rate: During conditioning, display the drift and add sample without starting the titration. Select a value below the highest drift as dosing rate. Stir faster.
	<i>The amount of methanol in the working medium is too low.</i>	<ul style="list-style-type: none"> Replace the working medium. Reduce the amount of solubility promoter, if working with solvent mixtures, see technical literature.
	<i>The electrode may be covered.</i>	Wipe off the electrode with ethanol or a suitable solvent.
The solution becomes darker after each titration.		Replace the working medium.
	<i>The electrode may be covered.</i>	Wipe off the electrode with ethanol or a suitable solvent.
	<i>The electrode has a short circuit.</i>	<ol style="list-style-type: none"> Check the Pt wires. Activate the electrode check.
The endpoint is reached too quickly.	<i>The dosing rate outside the control range is too high.</i>	Select the user-defined titration rate and reduce the dosing rate (see manual/help of the software used).
The titration times with volumetric titration are constantly longer.	<i>The buffer capacity of the solvent may be exhausted with two-component reagents.</i>	Replace the working medium.

7 Appendix

7.1 Remote interface

The 6.2148.010 Remote Box allows devices to be controlled which cannot be connected directly to the MSB interface of the Titrando.




Figure 24 Connectors of the Remote Box

1 Cable

For connecting to the Titrando.

2 MSB connector

Metrohm Serial Bus. For connecting external dosing devices or stirrers.

3 Remote connector

For connecting instruments with a remote interface.

7.1.1 Pin assignment of the remote interface





Figure 25 Pin assignment of remote socket and remote plug

The above figure of the pin assignment applies for all Metrohm instruments with 25-pin D-Sub remote connector.

Inputs




approx. 50 k Ω Pull-up

 $t_p > 20 \text{ ms}$

active = low, inactive = high

Outputs



Open Collector

 $t_p > 200 \text{ ms}$

active = low, inactive = high

$$I_C = 20 \text{ mA}, V_{CEQ} = 40 \text{ V}$$

+5 V: maximum load = 20 mA

The following tables offer information concerning the assignment of the individual pins and their function:

Table 6 Inputs and outputs of the remote interface

Assignment	Pin No.	Function*
Input 0	21	Start
Input 1	9	Stop
Input 2	22	
Input 3	10	Quit
Input 4	23	–
Input 5	11	
Input 6	24	
Input 7	12	
Output 0	5	Ready
Output 1	18	Conditioning OK
Output 2	4	Determination
Output 3	17	EOD
Output 4	3	
Output 5	16	Error
Output 6	1	
Output 7	2	Warning

Assignment	Pin No.	Function*
Output 8	6	
Output 9	7	
Output 10	8	
Output 11	13	
Output 12	19	
Output 13	20	
0 volts / GND	14	
+5 volts	15	
0 volts / GND	25	

* Signal activated only for operation with Touch Control.

Table 7 Explanation of the individual functions

Function	Explanation
Start	The current method is started at the time of activation. $t_{\text{pulse}} > 100 \text{ ms}$
Stop	The current method is canceled (Stop) at the time of activation. $t_{\text{pulse}} > 100 \text{ ms}$
Quit	The current command in the determination run will be canceled at the time of activation. $t_{\text{pulse}} > 100 \text{ ms}$
Ready	The instrument is ready to receive a start signal.
Conditioning OK	The line is set when Conditioning with SET titration and KFT titration is at OK. The line remains set until the determination is started with [START] .
Determination	The instrument performs a data-generating determination.
EOD	End of Determination. Pulse ($t_{\text{pulse}} = 200 \text{ ms}$) after a determination or after a buffer/standard solution during calibration using a Sample Processor.
Error	The line is set for error message display.

Function	Explanation
Warning	The line is set for warning message display.

8 Technical specifications

8.1 Measuring interface

The measuring cycle is 100 ms for all measuring modes.

8.1.1 Generator electrode

One measuring input (**Gen.**) for a generator electrode.

I_{max} 400 mA
Continuous current and pulsed

8.1.2 Indicator electrode

One measuring input (**Ind.**) for an indicator electrode.

Measuring mode Determination with adjustable polarization current
I_{pol}
AC 5, 10, 20 and 30 μ A
DC -125 - +125 μ A

8.1.3 Temperature

A measuring input (**Temp.**) for temperature sensors of the Pt1000 or NTC type with automatic temperature compensation.

R (25 °C) and B value can be configured for NTC sensors.

Measuring range
Pt1000 -150 - +250 °C
NTC -5 - +250 °C
(R (25 °C) = 30,000 Ω and B (25/50) = 4,100 K)

Resolution
Pt1000 0.1 °C
NTC 0.1 °C

Measuring accuracy
Pt1000 ± 0.2 °C
(Applies for measuring range -20 - +150 °C; ± 1 digit; without sensor error, under reference conditions)
NTC ± 0.6 °C
(Applies for measuring range +10 - +40 °C; ± 1 digit; without sensor error, under reference conditions)

8.5 Ambient temperature

Nominal function range	+5 - +45 °C
Storage	-20 - +60 °C
Transport	-40 - +60 °C

8.6 Interfaces

USB connectors

USB ports 2 USB downstream ports (type A sockets), 500 mA each, for connecting peripheral devices such as printers, keyboards, barcode readers or RS-232/USB boxes (Metrohm order no. 6.2148.020).

Connector "Controller"

<i>Controller port</i>	USB upstream port with auxiliary power supply (Mini DIN socket) for connecting Touch Control or computer for controlling the 852 Titrando.
<i>Touch Control</i>	With integrated Touch Control cable.
<i>Computer</i>	With 6.2151.000 cable.

MSB connectors (Metrohm Serial Bus)

<i>Dosing device</i>	
<i>Stirrer</i>	<p>Connection for a maximum of 4 stirrers.</p> <p>Stirrer control: switching on/off manually or coordinated with the titration sequence.</p> <p>Speed in 15 steps and shift direction can be selected.</p>
<i>Remote Box</i>	<p>Connection for a maximum of four Remote Boxes. Remote Boxes can be used to actuate and monitor external devices.</p>

9 Warranty (guarantee)

Metrohm guarantees that the deliveries and services it provides are free of defects in materials, design or manufacturing.

The general warranty period is 36 months (exclusions below) from the date of delivery or 18 months in the event of continuous operation. The warranty remains valid on the condition that the servicing is provided by a service organization authorized by Metrohm at defined intervals and with a defined scope.

The warranty period for anion suppressors of the type "MSM" is 120 months from the date of delivery or 60 months in the case of continuous operation.

The warranty period for IC separation columns is 90 days after start-up.

For third-party components that are recognizable as such, the manufacturer's warranty regulations apply.

Consumables and materials with limited storage life and glass breakage in the case of electrodes or other glass parts are excluded from the warranty.

Warranty claims cannot be asserted if the ordering party has failed to meet its payment obligations according to schedule.

During the warranty period, Metrohm undertakes either to replace free of charge or to credit the purchaser for any modules or components that can be shown to be faulty. Any transport or customs fees that may apply are the ordering party's responsibility.

The precondition for this is that the ordering party has to specify the article number, the article designation, an adequate error description, the delivery date and (if applicable) the serial number or chip data in the Support Tracker. Metrohm then decides whether a replacement or a credit note is to be issued or whether the faulty part has to be returned using the Return Material Authorization (RMA). If a replacement or credit note is issued, the ordering party undertakes to store the faulty part for at least 24 months in accordance with the current storage directives (in compliance with ESD guidelines) and to hold it in readiness for onsite inspection or for return shipment to Metrohm. Metrohm reserves the right to invoice the ordering party for these articles, including retroactively, in the event of noncompliance with these preconditions.

The same warranty periods that are specified for a corresponding new part apply to parts that are replaced or repaired within the above-mentioned warranty periods. However, replacement or repair of a part does not extend the warranty period of the entire system.

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty.

Metrohm also offers a 120 month spare parts availability guarantee and a 60 month PC software support warranty, calculated from the date on which the product is withdrawn from the market. The content of this warranty is the ability of the customer to obtain functioning spare parts or appropriate software support at market prices during the time of the warranty period.

If Metrohm AG is unable to meet this obligation due to circumstances beyond the control of Metrohm AG, then the ordering party shall be offered alternative solutions at preferential conditions.

10 Accessories



NOTE

Subject to change without notice.




10.1 Scope of delivery




10.1.1 2.852.0050





Qty.	Order no.	Description
1	1.801.0010	Magnetic stirrer
1	1.852.0010	Titrande
1	6.0338.100	Double Pt-wire electrode for volumetry Double Pt-wire electrode, Metrohm plug-in head G, for volumetric Karl Fischer titration and 684/737 Coulometer (cell with diaphragm). Shaft material: Glass Measuring range: -2000 - 2000 Measuring unit: mV Temperature range, long-term (°C): -20 - 70 Temperature range, short-term (°C): -20 - 70 SGJ sleeve: No Shaft diameter top (mm): 8 Shaft diameter bottom (mm): 8 Min. immersion depth (mm): 5 Electrode plug-in head: Metrohm plug-in head G Indicator electrode type: Pt Indicator electrode shape: Wire








Qty.	Order no.	Description	
1	6.0341.100	Double Pt wire electrode for coulometry	
		Double platinum wire electrode (0.8 x 4 mm), standard ground-joint, Metrohm plug-in head G, for Karl Fischer titrations.	
		Shaft material:	
		Glass	
		Measuring range:	
		–2000 - 2000	
		Measuring unit:	
		mV	
		Temperature range, long-term (°C):	
		–20 - 70	
		Temperature range, short-term (°C):	
		–20 - 70	
		SGJ sleeve:	
		Standard ground-joint 14/15	
		Shaft diameter top (mm):	12
		Shaft diameter bottom (mm):	8.75
		Min. immersion depth (mm):	10
		Electrode plug-in head:	Metrohm plug-in head G
		Indicator electrode type:	Pt
		Indicator electrode shape:	Wire
1	6.0344.100	Generator electrode with diaphragm	
		Generator electrode for Karl Fischer titrations, with diaphragm, standard ground-joint 29/22.	
		Shaft material:	
		Glass	
		Temperature range, long-term (°C):	
		–20 - 70	
		Temperature range, short-term (°C):	
		–20 - 70	
		SGJ sleeve:	
		Standard ground-joint 29/22	
		Shaft diameter top (mm):	24
		Shaft diameter bottom (mm):	24
		Min. immersion depth (mm):	15
		Electrode plug-in head:	Metrohm plug-in head G
		Indicator electrode type:	Cathode = Pt grid
		Indicator electrode type (reference):	Anode = Pt grid
1	6.1244.040	Set of seals	



Qty.	Order no.	Description	
1	6.1403.030	KF adsorber tube for coulometer cell	
		For KF coulometers.	
		Height (mm): 93	
		SGJ size: B-17/19	
1	6.1403.040	KF Adsorber tube	
		Used with KF equipment. Including cover and O-ring.	
		Height (mm): 113	
		Outer diameter (mm): 8	
1	6.1414.030	KF titration vessel lid	
		For volumetric KF titrations.	
		Material: PTFE	
		Material remark: Insert	

Qty.	Order no.	Description	
2	6.1448.010	Septum 12 mm, 5 pieces	
		Set of 5 items.	
		Material: Silicone rubber	
		Height (mm): 2	
		Outer diameter (mm): 12	
2	6.1448.020	Septum 16 mm, 5 pieces	
		Set of 5 pcs.	
		Material: Silicone rubber	
		Height (mm): 2	
		Outer diameter (mm): 16	
1	6.1464.320	KF titration vessel / 80-250 mL / coulometric	
		Material: Clear glass	
		Height (mm): 114	
		Volume (mL): 80 - 250	
2	6.1903.020	Stirring bar / 16 mm	
		Stirring Bar with magnetic core, PTFE coated, length 16 mm	
		Material: PTFE	
		Outer diameter (mm): 4	
		Length (mm): 16	




Qty.	Order no.	Description	
2	6.2016.030	Support rod / 250 mm Material: Stainless steel 18/9 Outer diameter (mm): 10 Length (mm): 250	
2	6.2043.005	Holding clip for bottles Holding spring for reagent bottles in exchange units.	
1	6.2047.020	Titration vessel holder Titration vessel holder for coulometric cells. Used with 756, 831 KF Coulometer	

Qty.	Order no.	Description	
1	6.2412.000	Weighing boats	
		Weighing boat made of glass for the 6.1414.030 titration vessel lid or the 6.1455.31X, 6.1464.32X, 6.1465.320 titration vessel. For KF water determinations. Including protective tube. Material: Glass Length (mm): 120	
1	6.2701.040	Screw cap	
		PBT. With orifice for 6.1448.020 Septum. Used with 6.5405.000 Cell without diaphragm and 6.1464.32X and 6.1465.320 Titration vessels. Used with 684, 737, 756 KF Coulometers Material: PBTP Height (mm): 20 Outer diameter (mm): 24 Inner diameter (mm): 18	
3	6.2713.000	Sleeve	
		For the fat-free ground-joint connection with the standard ground-joint SGJ 14. Material: PTFE foil	
1	6.2713.010	Sleeve	
		Sleeve for the fat-free ground-joint connection with the standard ground-joint SGJ 29. Material: PTFE foil	





Qty.	Order no.	Description
1	6.2738.000	Funnel
		For the KF coulometer.
		Material: PP
		Outer diameter (mm): 50
		Inner diameter (mm): 6
		
1	6.2811.000	Molecular sieve
		Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.
		
2	6.2816.030	Needle with Luer connector
		Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.
		Outer diameter (mm): 0.8
		Length (mm): 80
		




Qty.	Order no.	Description	
1	6.0338.100	Double Pt-wire electrode for volumetry	
		Double Pt-wire electrode, Metrohm plug-in head G, for volumetric Karl Fischer titration and 684/737 Coulometer (cell with diaphragm).	
		Shaft material:	
		Glass	
		Measuring range:	
		–2000 - 2000	
		Measuring unit:	
		mV	
		Temperature range, long-term (°C):	
		–20 - 70	
		Temperature range, short-term (°C):	
		–20 - 70	
		SGJ sleeve:	
		No	
		Shaft diameter top (mm):	
		8	
		Shaft diameter bottom (mm):	
		8	
		Min. immersion depth (mm):	
		5	
		Electrode plug-in head:	
		Metrohm plug-in head G	
		Indicator electrode type:	
		Pt	
		Indicator electrode shape:	
		Wire	
1	6.0341.100	Double Pt wire electrode for coulometry	
		Double platinum wire electrode (0.8 x 4 mm), standard ground-joint, Metrohm plug-in head G, for Karl Fischer titrations.	
		Shaft material:	
		Glass	
		Measuring range:	
		–2000 - 2000	
		Measuring unit:	
		mV	
		Temperature range, long-term (°C):	
		–20 - 70	
		Temperature range, short-term (°C):	
		–20 - 70	
		SGJ sleeve:	
		Standard ground-joint 14/15	
		Shaft diameter top (mm):	
		12	
		Shaft diameter bottom (mm):	
		8.75	
		Min. immersion depth (mm):	
		10	
		Electrode plug-in head:	
		Metrohm plug-in head G	
		Indicator electrode type:	
		Pt	
		Indicator electrode shape:	
		Wire	

Qty.	Order no.	Description	
1	6.1403.040	KF Adsorber tube	
		Used with KF equipment. Including cover and O-ring.	
		Height (mm): 113 Outer diameter (mm): 8	
1	6.1414.030	KF titration vessel lid	
		For volumetric KF titrations.	
		Material: PTFE Material remark: Insert	
1	6.1415.220	Titration vessel / 20-90 mL	
		Material: Clear glass	
		Height (mm): 80 Outer diameter (mm): 78 Volume (mL): 20 - 90	
1	6.1415.250	Titration vessel / 50-150 mL	
		Material: Clear glass	
		Height (mm): 80 Outer diameter (mm): 78 Volume (mL): 50 - 150	

Qty.	Order no.	Description
1	6.1464.320	KF titration vessel / 80-250 mL / coulometric
	Material:	Clear glass
	Height (mm):	114
	Volume (mL):	80 - 250
		
2	6.1903.020	Stirring bar / 16 mm
	Stirring Bar with magnetic core, PTFE coated, length 16 mm	
	Material:	PTFE
	Outer diameter (mm):	4
	Length (mm):	16
		
3	6.1903.030	Stirring bar / 25 mm
	Stirring bar with magnetic core, PTFE covering, length 25 mm.	
	Material:	PTFE
	Outer diameter (mm):	5
	Length (mm):	25
		




Qty.	Order no.	Description	
1	6.2713.010	Sleeve Sleeve for the fat-free ground-joint connection with the standard ground-joint SGJ 29. Material: PTFE foil	
1	6.2713.020	Sleeve For the fat-free ground-joint connection with the standard ground-joint SGJ 19. Material: PTFE foil	
2	6.2730.010	Screw nipple Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations. Material: PP Inner diameter (mm): 11.5 Length (mm): 23.5	
1	6.2730.020	Septum stopper With E.301.0041 O-ring. For 6.1448.010 Septum and 6.1414.030 Titration vessel lid. Used for KF applications. Outer diameter (mm): 18 Length (mm): 30	

Qty.	Order no.	Description
2	6.2816.030	Needle with Luer connector Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration. Outer diameter (mm): 0.8 Length (mm): 80
		
2	6.2816.090	Syringe 5 mL 5 mL syringe made of PP, with Luer connector. Material: PP
		
1	6.6064.010	USB Memory Stick for 900 Touch Control
		
1	6.2122.0x0	Power supply cable with IEC 60320 line socket, type C13 Cable plug according to customer requirements. Switzerland: SEV 1011, Type 12 6.2122.020 Germany, ...: CEE 7, Type VII 6.2122.040

Qty.	Order no.	Description
	USA, ...:	NEMA 5-15, Type 498 6.2122.070
	Length:	1.5 m
1	8.852.8003EN	852 Titrande Manual




10.1.3 2.852.0150





Qty.	Order no.	Description
1	1.852.0010	Titrande
1	6.0338.100	Double Pt-wire electrode for volumetry
<p>Double Pt-wire electrode, Metrohm plug-in head G, for volumetric Karl Fischer titration and 684/737 Coulometer (cell with diaphragm).</p> <p>Shaft material: Glass</p> <p>Measuring range: -2000 - 2000</p> <p>Measuring unit: mV</p> <p>Temperature range, long-term (°C): -20 - 70</p> <p>Temperature range, short-term (°C): -20 - 70</p> <p>SGJ sleeve: No</p> <p>Shaft diameter top (mm): 8</p> <p>Shaft diameter bottom (mm): 8</p> <p>Min. immersion depth (mm): 5</p> <p>Electrode plug-in head: Metrohm plug-in head G</p> <p>Indicator electrode type: Pt</p> <p>Indicator electrode shape: Wire</p>		






1	6.0341.100	Double Pt wire electrode for coulometry	
	Double platinum wire electrode (0.8 x 4 mm), standard ground-joint, Metrohm plug-in head G, for Karl Fischer titrations.		
	Shaft material:	Glass	
	Measuring range:	-2000 - 2000	
	Measuring unit:	mV	
	Temperature range, long-term (°C):	-20 - 70	
	Temperature range, short-term (°C):	-20 - 70	
	SGJ sleeve:	Standard ground-joint 14/15	
	Shaft diameter top (mm):	12	
	Shaft diameter bottom (mm):	8.75	
	Min. immersion depth (mm):	10	
	Electrode plug-in head:	Metrohm plug-in head G	



Qty.	Order no.	Description		
		Indicator electrode type:	Pt	
		Indicator electrode shape:	Wire	
1	6.0345.100	Generator electrode without diaphragm		
		Generator electrode for Karl Fischer titrations, without diaphragm, SGJ 29/22		
		Shaft material:		Glass
		Temperature range, long-term (°C):		-20 - 70
		Temperature range, short-term (°C):		-20 - 70
		SGJ sleeve:		SGJ 29/22
		Shaft diameter top (mm):		24
		Shaft diameter bottom (mm):		24
		Min. immersion depth (mm):		15
		Electrode plug-in head:		Metrohm plug-in head G
		Indicator electrode type:		Cathode = Pt sheet
		Indicator electrode type (reference):		Anode = Pt grid
1	6.1244.040	Set of seals		
1	6.1403.030	KF adsorber tube for coulometer cell		
		For KF coulometers.		
		Height (mm):		93
		SGJ size:		B-17/19



Qty.	Order no.	Description	
1	6.1437.000	SGJ stopper SGJ 14	
		Material:	Glass
		Height (mm):	34
		SGJ size:	B-14/15
1	6.1446.060	Stopper / B-14/15 / M10	
		Material:	PP
		Height (mm):	22
		SGJ size:	B-14/15
2	6.1448.010	Septum 12 mm, 5 pieces	
		Set of 5 items.	
		Material:	Silicone rubber
		Height (mm):	2
		Outer diameter (mm):	12
2	6.1448.020	Septum 16 mm, 5 pieces	
		Set of 5 pcs.	
		Material:	Silicone rubber
		Height (mm):	2
		Outer diameter (mm):	16



Qty.	Order no.	Description	
1	6.1464.320	KF titration vessel / 80-250 mL / coulometric	
		Material:	
		Height (mm):	
		Volume (mL):	
2	6.1903.020	Stirring bar / 16 mm	
		Stirring Bar with magnetic core, PTFE coated, length 16 mm	
		Material:	
		Outer diameter (mm):	
		Length (mm):	
3	6.1903.030	Stirring bar / 25 mm	
		Stirring bar with magnetic core, PTFE covering, length 25 mm.	
		Material:	
		Outer diameter (mm):	
		Length (mm):	

Qty.	Order no.	Description	
1	6.2001.090	Base plate on both sides Double stand for attaching 2 x 801, 803, 804 to a Titrande	
2	6.2013.010	Clamping ring For support rods with a diameter of 10 mm. Material: Metal Width (mm): 20 Height (mm): 16	
2	6.2016.030	Support rod / 250 mm Material: Stainless steel 18/9 Outer diameter (mm): 10 Length (mm): 250	
2	6.2043.005	Holding clip for bottles Holding spring for reagent bottles in exchange units.	


Qty.	Order no.	Description	
1	6.2104.120	Electrode cable / 1 m / H With plug H. For connecting generator electrode - KF Coulometer. Length (m): 1	
1	6.2412.000	Weighing boats Weighing boat made of glass for the 6.1414.030 titration vessel lid or the 6.1455.31X, 6.1464.32X, 6.1465.320 titration vessel. For KF water determinations. Including protective tube. Material: Glass Length (mm): 120	
1	6.2701.040	Screw cap PBT. With orifice for 6.1448.020 Septum. Used with 6.5405.000 Cell without diaphragm and 6.1464.32X and 6.1465.320 Titration vessels. Used with 684, 737, 756 KF Coulometers Material: PBTP Height (mm): 20 Outer diameter (mm): 24 Inner diameter (mm): 18	
3	6.2713.000	Sleeve For the fat-free ground-joint connection with the standard ground-joint SGJ 14. Material: PTFE foil	




Qty.	Order no.	Description
4	6.2730.030	Stopper
		Complete with nipple and E.301.0043 O-ring. For 6.1414.030 titration vessel lid. For KF instruments, VA measurement stands, VA Computrace and IC Dilution.
	Length (mm):	38
		
1	6.2738.000	Funnel
		For the KF coulometer.
	Material:	PP
	Outer diameter (mm):	50
	Inner diameter (mm):	6
		
1	6.2811.000	Molecular sieve
		Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.
		

10.1.4 2.852.0160





Qty.	Order no.	Description
1	1.852.0010	Titrande
1	1.900.0010	Touch Control
1	6.0338.100	Double Pt-wire electrode for volumetry Double Pt-wire electrode, Metrohm plug-in head G, for volumetric Karl Fischer titration and 684/737 Coulometer (cell with diaphragm). Shaft material: Glass Measuring range: -2000 - 2000 Measuring unit: mV Temperature range, long-term (°C): -20 - 70 Temperature range, short-term (°C): -20 - 70 SGJ sleeve: No Shaft diameter top (mm): 8 Shaft diameter bottom (mm): 8 Min. immersion depth (mm): 5 Electrode plug-in head: Metrohm plug-in head G Indicator electrode type: Pt Indicator electrode shape: Wire
1	6.0341.100	Double Pt wire electrode for coulometry Double platinum wire electrode (0.8 x 4 mm), standard ground-joint, Metrohm plug-in head G, for Karl Fischer titrations. Shaft material: Glass Measuring range: -2000 - 2000 Measuring unit: mV Temperature range, long-term (°C): -20 - 70 Temperature range, short-term (°C): -20 - 70 SGJ sleeve: Standard ground-joint 14/15 Shaft diameter top (mm): 12 Shaft diameter bottom (mm): 8.75 Min. immersion depth (mm): 10 Electrode plug-in head: Metrohm plug-in head G Indicator electrode type: Pt Indicator electrode shape: Wire






Qty.	Order no.	Description	
1	6.1403.040	KF Adsorber tube	
		Used with KF equipment. Including cover and O-ring.	
		Height (mm): 113 Outer diameter (mm): 8	
1	6.1414.030	KF titration vessel lid	
		For volumetric KF titrations.	
		Material: PTFE Material remark: Insert	
1	6.1415.220	Titration vessel / 20-90 mL	
		Material: Clear glass	
		Height (mm): 80 Outer diameter (mm): 78 Volume (mL): 20 - 90	
1	6.1415.250	Titration vessel / 50-150 mL	
		Material: Clear glass	
		Height (mm): 80 Outer diameter (mm): 78 Volume (mL): 50 - 150	

Qty.	Order no.	Description
1	6.1464.320	KF titration vessel / 80-250 mL / coulometric
	Material:	Clear glass
	Height (mm):	114
	Volume (mL):	80 - 250
		
2	6.1903.020	Stirring bar / 16 mm
	Stirring Bar with magnetic core, PTFE coated, length 16 mm	
	Material:	PTFE
	Outer diameter (mm):	4
	Length (mm):	16
		
3	6.1903.030	Stirring bar / 25 mm
	Stirring bar with magnetic core, PTFE covering, length 25 mm.	
	Material:	PTFE
	Outer diameter (mm):	5
	Length (mm):	25
		

Qty.	Order no.	Description	
1	6.2047.020	Titration vessel holder Titration vessel holder for coulometric cells.Used with 756, 831 KF Coulometer	
2	6.2103.130	Adapter red 2 mm plug / 4 mm socket For connecting plug B (4 mm) to 2 mm socket.	
2	6.2103.140	Adapter black 2 mm plug / B socket 4 mm For connecting plug B (4 mm) to 2 mm socket.	
2	6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F). Length (m): 1	


Qty.	Order no.	Description	
1	6.2713.010	Sleeve Sleeve for the fat-free ground-joint connection with the standard ground-joint SGJ 29. Material: PTFE foil	
1	6.2713.020	Sleeve For the fat-free ground-joint connection with the standard ground-joint SGJ 19. Material: PTFE foil	
2	6.2730.010	Screw nipple Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations. Material: PP Inner diameter (mm): 11.5 Length (mm): 23.5	
1	6.2730.020	Septum stopper With E.301.0041 O-ring. For 6.1448.010 Septum and 6.1414.030 Titration vessel lid. Used for KF applications. Outer diameter (mm): 18 Length (mm): 30	

Qty.	Order no.	Description
2	6.2816.030	Needle with Luer connector
		Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.
		Outer diameter (mm): 0.8
		Length (mm): 80
		
2	6.2816.090	Syringe 5 mL
		5 mL syringe made of PP, with Luer connector.
		Material: PP
		
1	6.6064.010	USB Memory Stick for 900 Touch Control
		
1	6.2122.0x0	Power supply cable with IEC 60320 line socket, type C13
		Cable plug according to customer requirements.
		Switzerland: SEV 1011, Type 12 6.2122.020
		Germany, ...: CEE 7, Type VII 6.2122.040



Qty.	Order no.	Description
	USA, ...:	NEMA 5-15, Type 498 6.2122.070
	Length:	1.5 m
1	8.852.8003EN	852 Titrande Manual

10.2 Optional accessories

Order no.	Description
2.136.0100	Polytron PT 1300 D
Polytron PT 1300 D - Metrohm Version Homogenizer with digital display that can be actuated directly via tiamo™ or Touch Control.	
2.800.0010	800 Dosino
Drive with write/read hardware for intelligent Dosing Units. With fixed cable (length 150 cm).	
2.801.0010	801 Stirrer
Magnetic stirrer without stand for supplementing the Titrino plus, Dosimat plus, Titrandos, Sample Processors, 805 Dosimats and 780/781 pH Meters. With permanently attached cable for MSB (Metrohm Serial Bus).	
2.803.0010	803 TI Stand with stirrer and pump
Compact titration stand with magnetic stirrer and built-in pump for solvent addition and aspiration of the titration vessel contents. Suitable for Karl Fischer titrations with Titrandos. MSB connection. Including comprehensive accessories, solvent bottle and waste bottle.	
Dimensions in mm (W/H/D):	106/ 101/ 220
Dimensions remark:	Height with support rod: 412 mm

Order no.	Description	
2.860.0010	860 KF Thermoprep	
<p>The 860 KF Thermoprep is designed for thermal sample preparation in Karl Fischer titration. Numerous substances cannot be analysed by direct Karl Fischer titration as they are not soluble, react with the Karl Fischer reagent or release their water only very slowly or not until heated to high temperatures. Using tightly sealed sample vials, the sample is inserted in the oven. The samples can be subsequently analysed by any volumetric or coulometric KF titrator.</p> <p>Dimensions in mm (W/H/D): 440/500/230</p>		
6.2148.010	Remote Box MSB	
<p>Additional remote interface for the connection of devices that can be controlled via remote lines. With permanently attached cable.</p>		
6.2148.020	RS-232/USB box	
<p>Extends the Titrand, USB Sample Processors, Dosing Interface, Titrosampler by two RS-232 interfaces for connecting e.g. balances, PC, LIMS. Requires 6.2151.030 cable.</p>		
6.3032.120	Dosing Unit 2 mL	
<p>Dosing unit with integrated data chip with 2 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 glass thread. FEP tubing connection, anti-diffusion buret tip.</p> <p>Volume (mL): 2</p>		
6.3032.150	Dosing Unit 5 mL	
<p>Dosing unit with integrated data chip with 5 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.</p> <p>Volume (mL): 5</p>		
6.3032.210	Dosing Unit 10 mL	
<p>Dosing unit with integrated data chip with 10 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.</p> <p>Volume (mL): 10</p>		

Order no.	Description	
6.3032.220	Dosing Unit 20 mL	
<p>Dosing unit with integrated data chip with 20 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.</p> <p>Volume (mL): 20</p>		
6.3032.250	Dosing Unit 50 mL	
<p>Dosing unit with integrated data chip with 50 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, anti-diffusion buret tip.</p> <p>Volume (mL): 50</p>		
6.5406.000	Coulometric KF equipment, including generator electrode with diaphragm	
<p>Complete accessory set for KF coulometer.</p>		
6.5406.010	Coulometric KF equipment, including generator electrode without diaphragm	
<p>Complete accessory set for KF coulometer.</p>		
6.5617.000	Equipment for reagent replacement	
<p>Equipment for automatic reagent replacement with KF coulometer and connected Dosino.</p>		
2.874.0010	874 USB Oven Sample Processor	
<p>The 874 Oven Sample Processor is used for thermal sample preparation in Karl Fischer titration. The oven method is suitable for samples that only release their moisture at higher temperature, for sparingly soluble samples or samples that undergo side-reactions with the KF reagent.</p> <p>Dimensions in mm (W/H/D): 280/550/490</p>		

Order no.	Description	
2.874.0020	874 USB Oven Sample Processor customized	
<p>The 874 Oven Sample Processor is used for thermal sample preparation in Karl Fischer titration. The oven method is suitable for samples that only release their moisture at higher temperature, for sparingly soluble samples or samples that undergo side-reactions with the KF reagent.</p> <p>Dimensions in mm (W/H/D): 280/550/490</p>		
2.900.0010	900 Touch Control	
<p>Operating unit for the Titrandos, USB Sample Processors, 856 Conductivity Module, 867 pH Module and 846 Dosing Interface. Touch-sensitive, high-resolution color display, simple and intuitive operation, thanks to favorites for direct method access. With integrated Ethernet interface for direct connection to the Internet and USB interface for connecting USB printers or a USB memory stick.</p>		
6.6056.231	tiamo™ 2.3 Light CD: 1 license	
<p>tiamo™ 2.3 light computer program for controlling a titration system. Up to two instruments can be connected Graphical method editor with numerous templates Layout Manager for customizing the user interface Professional database with reevaluation High-performance report generator No parallel titration, no data export 1 license Dialog languages: German, English</p>		
6.6056.232	tiamo™ 2.3 Full CD: 1 license	
<p>computer program for controlling complex titration systems.</p>		
6.6056.233	tiamo™ 2.3 Multi CD: 3 licenses	
<p>Client/server-compatible computer program for controlling complex titration systems.</p>		

Numbers/Symbols

Reagent replacement 46, 51

Remote	
Interface	59
Pin assignment	59
Remote Box	
Connect	18

S

Safety instructions	5
Sample	
Liquid	44, 50
Solid	45, 50
Sample changer	
With oven module	33
Sample size	
Size	44, 48
Sensor	
Connect	38
Serial number	9
Service	6
SET	3
Stirrer	
Connect	17
Stirring direction	24
System test	12

T

Temperature sensor	
Connect	38

Titrand system	1
Titration cell	
Equip	27
Fill	28
Titration mode	3
BRC	3
KFC	3
KFT	3
MET	3
SET	3
Titration stand	
Connect	17
Titration vessel	
Equipping	24
Touch Control	
Connect	12

U

Update	
Device software	2
USB	
Connector	9
USB connector	2
USB hub	
Connect	20

V

Validation	54
------------------	----

Volumetric KF titration	
Application	48
Endpoint determination, biam-	
perometric (Upol)	47
Endpoint determination, bivol-	
tametric (Ipol)	47
One-component reagents ...	48
Principle	47
Sample size	48
Tips for daily practice	48
Two-component reagents .	48
Water standard	48
Working conditions	51
Volumetric titration	
KFT	4
MET	4
SET	4

W

Warranty	67
Water standard	
Certified	41, 48